

DELUXE CUBRO

USER GUIDE FOR DELUXE® FOLDING WALKING FRAME



MAX USER WEIGHT: 120KG

ORDERING CODE: 1800, 1801, 1802, 1803, 1804, 1805, 1806

CONTENTS

1. Intended Use	2
2. Product Description	2
3. Definitions	2
4. Warnings	3
5. Product Specifications	4
6. List of Components	4
7. Instructions for Use	5
8. Care and Maintenance	5
9. Limited Warranty	7

1. INTENDED USE

The Deluxe folding walking frame is designed for individuals who, due to a wide variety of reasons, require additional support while walking. It should be used as a tool to assist people who have mobility difficulties with standing, walking or sitting independently.

2. PRODUCT DESCRIPTION

Manufactured from high quality materials these lightweight aluminium walking frames are height adjustable and fold away for easy storage. The hoop design gives maximum strength, stability and long life.

3. DEFINITIONS

Symbols used in this user manual and their meanings.



WARNING! Failure to heed this warning may result in damage to the product or serious injury to the operator/user.



CAUTION! Beware of pinch points.



ATTENTION! Read and understand the instructions in the user guide before using the product.

4. WARNINGS



ATTENTION! This user manual must be read before using the walking frame. Severe injury or death may result if user instructions, maintenance instructions and product warnings are not followed.



WARNING! The decision to use this walking frame should always be guided by a qualified Health Care Professional, who can also advise and assist with the correct height adjustment and usage instruction.



CAUTION! Fingers or other body parts can be trapped between the walking frame's folding mechanism during operation. Be extremely careful when adjusting the walking frame. Observe the dimensions of gaps between the walking frame's elements to avoid trapping your fingers or other body parts.


LIMITATIONS FOR USE

- Maximum user weight is 120kg
- Instruction and practice are necessary for successful use of the walking frame.
- The walking frame is to be used indoors on flat surfaces.
- Use the walking frame only for the purpose it is intended for.
- DO NOT stand on the walking frame.
- DO NOT hang anything from the frame of your unit.
- Ensure legs, wheels, stops, skis are level to avoid instability problems.
- DO NOT use to pull-up on. This is not a transfer aid.

5. PRODUCT SPECIFICATIONS

SPECIFIC DIMENSIONS

Maximum user weight	120kg
Adjustable height	790mm - 970mm
Width between handles	440mm
Overall width	615mm
Inside width	560mm
Tip size	28mm
Ski size	25mm



6. LIST OF COMPONENTS

The Deluxe Folding Walking Frames are supplied in boxes and come folded, for easy transportation. The contents include one complete Deluxe Folding Walking Frame and one user manual.



OPTIONAL
Basket and Tray
Ordering Code: 1860

Description	4 rubber feet	Front wheels, rear stops	4 fixed wheels	Front castors, rear stops	Front wheels, rear skis	Front castors, rear wheels	Front castors, rear skis
Code	1800	1801	1802	1803	1804	1805	1806
Front							
Rear							
Frame Weight	2.5kg	3.3kg	4.1kg	3.6kg	3.3kg	4.2kg	3.5kg

7. INSTRUCTIONS FOR USE

TO UNFOLD

- Swing side frames outwards until you hear the spring button “click” into place.
- Check frame is locked and rigid before use.



WARNING! Ensure the frame is fully open before use.

TO FOLD


- Depress button (under black plastic cover) in centre of top folding tube.
- Swing side frames inwards.



CAUTION! Beware of pinch points

HEIGHT ADJUSTMENT

- Ensure left and right legs are located in the same holes and spring buttons protrude through adjustment holes in stems.



WARNING! Seek advice from a qualified health professional for the appropriate height level to suit your needs.

8. CARE AND MAINTENANCE

CLEANING

- Clean the walking frame regularly with a damp cloth
- Use a mild washing agent and warm water for excess dirt build up.
- Do not use harsh abrasives or bleach based agents which may cause damage or discolouration.

WHEELS

- Clean wheels, castors, stops and skis with a mild washing agent and warm water.
- Lubricate wheels with multi-purpose oil periodically or if wheels begin to squeak.

MAINTENANCE

- Ensure that every part of the walking frame is secure at all times.
- Check all screws for tightness. If not, tighten if necessary.
- Replace any broken, damaged or worn items immediately.

DISPOSAL

Products that can no longer be used are to be disposed of separately from household rubbish. This must be done according to the local and national regulations for environmental protection and raw material recycling. Please look at your local council's website for further information on the correct disposal

MAINTENANCE SCHEDULE

Weekly	Periodically	Yearly	Serviceable Components
<ul style="list-style-type: none"> • Check all wheels, nuts and bolts are secure 	<ul style="list-style-type: none"> • Check handgrips for wear • Check folding hinge bolts are secure • Check front and rear axle bolts are secure • Check frame for wear or cracks • Check wheels for damage or splitting 	<ul style="list-style-type: none"> • Take your walking frame to an authorised dealer for a full check, clean and service. 	<ul style="list-style-type: none"> • Hand grips • Wheels • Bearings • Nuts and bolts



WARNING! If faults found during any inspection remove from use immediately.



WARNING! Do not make any changes or modifications to this product without consultation from your Cubro agent.

For Any Servicing, Maintenance and Troubleshooting Queries,
Please Contact The Cubro Service Dept;

Ph. 0800 656 527

Ph. +64 7 578 7228

E: sales@cubro.co.nz

cubro.co.nz

Address: 149 Taurikura Dr, Tauriko, Tauranga, New Zealand

9. LIMITED WARRANTY

1. This warranty applies only to the original purchaser of CUBRO LTD's product (who must be a purchaser who is resident or carrying on business in New Zealand and who has purchased the product directly from CUBRO LTD, or from a CUBRO LTD authorized distributor or reseller in New Zealand) and it is not transferable to any other person or entity.
2. This warranty replaces all conditions, warranties or guarantees that might be implied by law in favour of any person, which are excluded to the fullest extent permissible by law. Where the Consumer Guarantees Act would apply but the purchaser acquires or holds themselves out as acquiring any product for business purposes, the guarantees provided under the Consumer Guarantees Act are excluded in relation to that product.
3. CUBRO LTD warrants the product to be free from defects in materials and workmanship from date of purchase for a period of one year (12 months). This warranty does not cover any damage, defect, expense or loss of any kind caused by accident, misuse, abuse, neglect, negligence, alteration or modification (which includes the use of unauthorised parts or attachments), improper service, repair by other than authorized personnel or any defects not related to materials or workmanship. Wear of components in normal operation and failures resulting therefrom are excluded from this warranty.
4. CUBRO LTD has no obligation to the purchaser and is not obligated to honour all or any part of this warranty unless the following procedure is followed by the purchaser:
 - a. Before making a warranty claim, the purchaser should ensure that the product is defective by following standard 'trouble shooting' procedures and be able to attribute the fault to a defect in materials and/or workmanship of CUBRO LTD;
 - b. If CUBRO LTD determines that the product is not defective, or that the warranty claim is otherwise invalid, CUBRO LTD shall charge the purchaser a handling and servicing fee as well as any freights costs.
 - c. Upon receiving notice from the purchaser of an alleged defect in a product (which defect should be reported to CUBRO LTD immediately), CUBRO LTD will issue a serialized return management authorization (RMA). The purchaser must then return the entire unit or remove, at the purchaser's cost, the defective component part(s) identified, pack the unit or the component part(s) in a manner to avoid shipping damage and to ship the unit or the component part(s) to either CUBRO LTD, or a service centre as specified by CUBRO LTD, within 30 days of the date of the serialized return authorization date.
 - d. If CUBRO LTD require additional information relating to the use of the product by the purchaser or any other relevant information, the purchaser will supply such information as soon as practically possible and in such a manner as CUBRO LTD may reasonably require.
 - e. If CUBRO LTD access the purchaser's claim, it will either repair or replace the product, or repair or replace the defect in workmanship, as it may determine in its sole discretion.
5. Proper selection of a specific product for a specific application and operating environment, and its compatibility with other equipment is the purchaser's responsibility. CUBRO LTD does not warrant the performance of its products or their suitability for a particular purpose.
6. CUBRO LTD shall not be liable under any heading (whether in contract, negligence, or otherwise) for any consequential, indirect or incidental loss or damages of any kind and the maximum aggregate liability of CUBRO LTD (under any one or more headings) shall be the purchase price of the product.