

Shower commode

The Etac® Clean shower commode is an ideal solution for those wanting a secure, safe and simple option for both users and carers.

It provides a comfortable fit for many individuals and is simple to manoeuvre.

With wide push handles and a comfortable seating area, the Etac Clean shower commode is packed full of features which provide additional support.

Comes with:



Pan holder



Pan with lid

Optional accessories:



Padded seat
2685



Safe working load		130kg
Width between arm supports		435mm
Seat height		550mm
Seat dimensions	W	480mm
	D	440mm
Overall dimensions	W	500mm
	D	565mm
	H	1070mm
Toilet bowl clearance	with pan holder	480mm
	without pan holder	495mm
Seat aperture	W	360mm
	D	215mm
Product weight		14.5kg

Features

- Compact overall dimensions while maintaining a large and comfortable seating area
- Upright operating handle that allows the chair to fit over most toilets, including wall-mounted
- Adjustable foot support has a soft, rounded shape which gives the feet extra comfort and protection
- Single tube frame has no welded joints to reduce the chance of rust
- Rear opening seat to facilitate easier hygiene access
- Powder-coated steel for extra durability
- 4x 125mm lockable castors

Code: 2680

Etac® Clean shower commode



A safe, simple and secure choice for users and carers



The Etac® Clean shower commode is an ideal solution for those wanting a secure, safe and simple option for both carers and users. Practical functionality and design combine to create a positive bathroom and showering experience.

The upright handle ensures that Clean is a perfect fit, even over wall-mounted toilets, and the chair is open at the back to facilitate access.

With wide push handles, the Etac® Clean is easy for carers to move around a shower room. It's also designed to be easily washed down and cleaned, with a no-weld design and enclosed joins and castors.

The commode comes complete with a pan with lid and pan holder.