

TRAVELLER

Transit wheelchair

CUBRO TRAVELLER

Ultra lightweight for easy transportation

The TRAVELLER Transit wheelchair is an ultra lightweight wheelchair for easy transportation. The TRAVELLER also features an ergonomic design for ease of use.

Features

- Quick release fold down backrest
- Swing-away, length adjustable footplates
- Full length height adjustable armrests
- Curb climber feature
- Adjustable footplate height
- 2x solid 12" (305mm) rear wheels and 2x 8" (200mm) front castors

Code: 3046



Specifications



Seat width

460mm



Overall width

700mm



Overall depth

1060mm



Overall height

930mm



Folded width

300mm



Seat depth

410mm



Seat height

510mm



Armrest height

210mm



Back height

450mm



Product weight

11.3kg



09 283 5930



www.phillipshealthcare.co.nz