

Heel Protector Boot

POSIMED®

As part of systematic decubitus therapy, the Posimed® Heel Protector Boot is a medical device constructed of breathable open cell Geltouch® and high resilience (HR) foam.

It is designed to offload pressure from the heel of a non-ambulatory individual by elevating and cushioning the heel, in turn reducing pressure, friction and shear from the surface of a bed or wheelchair, and subsequent pressure injuries.





Autoclavable according to RKT® guidelines

Suitable for the treatment of diabetic foot ulcers

For decubitus prophylaxis and therapy

Cover is easy to disinfect



Features

- Easy to disinfect by wiping
- Washable up to 95°C
- Suitable for tumble dryer
- Autoclavable

Code: 431