




# Vermeiren Jazz light manual wheelchair



CRASH TESTED



**Lightweight and highly portable**



**Self-propelled wheelchair**



**Ergonomic hand grips**

The Vermeiren Jazz manual wheelchair is a self propelled wheelchair, it has been crash tested and is made to last. Features include ergonomic handgrips, removable armrests and adjustable legrests. Easily folds down for efficient storage and travel.

Swing away footrests (both inward and outwards) and foldable footplates for ease of transfers.

## Models

Vermeiren Jazz light manual wheelchair 460mm wide

564-46

Vermeiren Jazz light manual wheelchair 500mm wide

564-50

## Codes

## SPECIFICATIONS

 1120mm L 935mm H Overall dimensions	 130kg Safe working load
 220mm Armrest height from seat	 420mm Backrest height from seat
 Crash tested Compliance	 min. 370mm max. 440mm Footrest adjustment
 935mm Max height	 16.6kg Product weight
 Belgium Made in	

## FEATURES

- Ergonomic handgrips
- Removable armrests
- Self-propelled folding wheelchair
- Lightweight and highly portable
- Swing away footrests (both inward and outwards) and foldable footplates for ease of transfers
- Adjustable and removable legrests
- 200mm solid PU front wheels and 610mm solid PU rear wheels
- ISO 7176-19 compliant



09 283 5930



[www.phillipshealthcare.co.nz](http://www.phillipshealthcare.co.nz)