

# Milano® Overbed table

The Milano® Overbed table features a one-piece, easy-to-clean top designed for superior infection control. Suitable for use over both beds and chairs, it is easily height adjustable and its rounded edge helps to prevent items from falling off while the user is eating or reading.



## Features

- One-piece, moulded ABS plastic top for superior infection control
- Soft rounded edge helps prevent items from slipping off
- Gas-assisted height adjustment with a left hand handle that allows the table to be easily adjusted with one hand
- Fits over most beds and chairs
- Easy clean and care
- Modern and robust
- 4x 50mm twin castors, 2x braking

**Code: 594**



Maximum table load	20kg
Table dimensions	L 900mm W 430mm
Height range (including castors)	725mm - 1050mm
Width inside base	750mm
Product weight	14.6kg