

# Mobilis® walking frames

Mobilis® walking frames are designed to assist individuals who have mobility difficulties with standing, walking or sitting independently. Their aluminium frames mean they are lightweight, easy to transport and to store. Their superior handbrakes are 'one touch' and easy to adjust. They all come standard with a padded seat and underseat storage bag.



## The Mobilis® Standard selection

### Mobilis® walking frame

Our standard Mobilis® walking frame is a lightweight frame with 150mm front castors and rear wheels making it suitable for indoor use.



## The Mobilis® Narrow range

### Mobilis® Narrow walking frame

This frame is shorter, narrower and lower than most standard walking frames meaning it is perfect for use in space restricted areas and for people of smaller stature. It has 150mm front castors and rear wheels making it suitable for indoor use.



Extra tall handles may be available on some of these ranges. Talk to us about these options today.



**NEW**

### Mobilis® Narrow Super Stroller-8 walking frame

This frame is shorter, narrower and lower than most standard walking frames meaning it is perfect for use in space restricted areas and for people of smaller stature. Its larger 200mm (8") front castors and rear wheels make it suitable for indoor and outdoor use.

## Get more from your Mobilis®

### Mobilis® Quad walking frame

Has the same features as the standard Mobilis® walking frame but with larger 200mm (8") front castors and rear wheels making it suitable for indoor and outdoor use.

Also available with optional extras including a folding tray and extra tall handles.



### Mobilis® Plus walking frame

The Mobilis® Plus has all the features of the Mobilis® Quad frame and more. Extra features include a cane holder that enables users to transport their cane on their walking frame and a height adjustable seat.

Available with optional extra of a folding tray.




All Mobilis® walking frames fold down easily for transport and storage



# CUBRO Mobilis® walking frame range

## Mobilis® walking frames at a glance

A quick and easy comparison of our Mobilis® range.



Model	Colours	Outdoor use	SWL	Product weight	Overall dimensions	Inside width	Seat size	Handle height
<b>The Mobilis® Narrow range</b>								
<b>Mobilis® Narrow walking frame</b> <ul style="list-style-type: none"> <li>Shorter, narrower and lower frame</li> <li>Ideal for space restricted areas</li> <li>Suits smaller statures</li> <li>150mm front castors and rear wheels</li> </ul>	Silver	No	130kg	6kg	W 500mm D 600mm	445mm	W 300mm D 300mm H 470mm	Min 735mm Max 850mm
<b>Mobilis® Narrow Super Stroller-8 walking frame</b> <ul style="list-style-type: none"> <li>Shorter, narrower and lower frame</li> <li>Ideal for space restricted areas</li> <li>Suits smaller statures</li> <li>200mm front castors and rear wheels</li> </ul>	Aqua	Yes	120kg	7.6kg	W 540mm D 670mm	390mm	W 305mm D 305mm H min. 450mm max. 500mm	Min 745mm Max 840mm
<b>The Mobilis® Standard selection</b>								
<b>Mobilis® walking frame</b> <ul style="list-style-type: none"> <li>Lightweight frame</li> <li>150mm front castors and rear wheels</li> </ul>	Blue Red Silver	No	120kg	6.5kg	W 610mm D 600mm	445mm	W 305mm D 305mm H 530mm	Min 820mm Max 930mm
<b>The Mobilis® Quad and Mobilis® Plus</b>								
<b>Mobilis® Quad walking frame</b> <ul style="list-style-type: none"> <li>200mm front castors and rear wheels</li> <li>Optional folding tray and extra tall handles</li> </ul>	Green	Yes	120kg	7.5kg	W 610mm D 660mm	445mm	W 350mm D 350mm H 530mm	Min 860mm Max 920mm
<b>Mobilis® Plus walking frame</b> <ul style="list-style-type: none"> <li>200mm front castors and rear wheels</li> <li>Cane holder</li> <li>Height adjustable seat</li> <li>Optional folding tray</li> </ul>	Red Silver Blue Flame Green Flame Purple	Yes	120kg	7.5kg	W 620mm D 660mm	435mm	W 350mm D 300mm H min. 520mm max. 580mm	Min 760mm Max 960mm