




.bock[®]///

USER GUIDE

INTENDED USE


The clamp on support side rail assists with bed mobility, enables the use of both hands to sit up on the side of the bed and then also offers a means to maintain balance while rising to a standing position. Support handles are easy to attach and remove, and are designed to fit over 20mm bed frames. The safe working load for this side rail is 150kg.


DIMENSIONS




DEFINITIONS

Symbols used in this user manual and their meanings:

 **WARNING!** Failure to heed this warning may result in damage to the product or serious injury to the operator/user.

 **CAUTION!** Beware of entrapment.

 **ATTENTION!** Read and understand the instructions in the user guide before using the product.

Symbols used on the support rail and their meanings:



- SAFE WORKING LOAD 150KG.
- WARNING Incompatible SIDE RAILS can create HAZARDS. Consult instructions for use to confirm compatibility.